




# Advance<sup>®</sup> 14LP

LOW PROFILE PTA BALLOON DILATATION CATHETER



True peripheral design – specifically for infrapopliteal PTA.

#### PERIPHERAL SHAFT DESIGN

Coiled wire shaft provides unmatched pushability and trackability while delivering flexibility and kink resistance.

#### LOW-PROFILE SYSTEM

Ultra-low crossing profile with 4 Fr sheath compatibility on all sizes, even after inflation.

#### LONG LENGTH BALLOONS

Treat infrapopliteal lesions with balloon lengths from 20 to 200 mm.

LABORATORY  
TESTING\*

\*Data on file.

ADVANCING ACCESS > TARGET > TREAT

LEG THERAPIES

True peripheral  
**design** makes the most  
challenging infrapopliteal  
procedures possible.





Unlike balloons designed for coronary applications and repurposed for peripheral use, the .014" wire compatible Advance 14LP has been designed specifically to combat infrapopliteal disease with our exclusive **coiled wire shaft**. Why use a coronary balloon for below-the-knee PTA when you can use one that was designed for it?



Advance® 14LP

LOW PROFILE PTA BALLOON DILATATION CATHETER

True peripheral design translates into superior performance.



Amphirion is a registered trademark of Invatec S.R.L.

Sleek is a registered trademark of Cordis Corporation.

Sterling is a registered trademark of Boston Scientific Scimed, Inc.

## PERIPHERAL SHAFT DESIGN

### Exclusive Coiled Wire Shaft

- Provides the trackability and pushability to reach even the most remote infrapopliteal lesions.
- Delivers flexibility and kink resistance, which enhance your ability to travel through the tortuous vasculature of the lower leg.
- Ensures stable wire position during sheath retraction.

## LOW PROFILE SYSTEM

### Ultra-Low Crossing Profile

- Hydrophilic coating on balloon and distal shaft, along with a smooth tip transition, improves ability to cross tight lesions and reduces resistance.


### 4 Fr Sheath Compatibility

- Maintains super-low profile after inflation, facilitating sheath pullback and reinsertion – plus 4 Fr sheath compatibility for all sizes.

## LONG LENGTH BALLOONS

### Balloon Sizes up to 200 mm

- Treat longer lesions with fewer inflations and deflations with balloon sizes from 20 to 200 mm long in 2, 2.5, 3 and 4 mm diameters.



Exclusive coiled wire shaft for unmatched performance



# Advance® 14LP


LOW PROFILE PTA BALLOON DILATATION CATHETER

Global Product Number	Order Number	Compatible Sheath Fr	Nominal Pressure atm	Rated Burst Pressure atm
<b>170 cm Catheter Length</b>				
G50313	PTAX4-14-170-2-2	4	8	12
G50314	PTAX4-14-170-2-4	4	8	12
G50315	PTAX4-14-170-2-6	4	8	12
G50316	PTAX4-14-170-2-8	4	8	12
G50317	PTAX4-14-170-2-12	4	8	12
G50318	PTAX4-14-170-2-16	4	8	12
G50319	PTAX4-14-170-2-20	4	8	12
G50320	PTAX4-14-170-2.5-2	4	8	12
G50321	PTAX4-14-170-2.5-4	4	8	12
G50322	PTAX4-14-170-2.5-6	4	8	12
G50323	PTAX4-14-170-2.5-8	4	8	12
G50324	PTAX4-14-170-2.5-12	4	8	12
G50325	PTAX4-14-170-2.5-16	4	8	12
G50326	PTAX4-14-170-2.5-20	4	8	12
G50327	PTAX4-14-170-3-2	4	8	12
G50328	PTAX4-14-170-3-4	4	8	12
G50329	PTAX4-14-170-3-6	4	8	12
G50330	PTAX4-14-170-3-8	4	8	12
G50331	PTAX4-14-170-3-12	4	8	12
G50332	PTAX4-14-170-3-16	4	8	12
G50333	PTAX4-14-170-3-20	4	8	12
G50334	PTAX4-14-170-4-2	4	8	12
G50335	PTAX4-14-170-4-4	4	8	12
G50336	PTAX4-14-170-4-6	4	8	12
G50337	PTAX4-14-170-4-8	4	8	12
G50338	PTAX4-14-170-4-12	4	8	12
G50339	PTAX4-14-170-4-16	4	8	12
G50340	PTAX4-14-170-4-20	4	8	12

The Advance 14LP PTA Balloon Dilatation Catheter has been designed for percutaneous transluminal angioplasty (PTA) of lesions in peripheral arteries including iliac, renal, popliteal, infrapopliteal, femoral and iliofemoral as well as obstructive lesions of native or synthetic arteriovenous dialysis fistulas.

Please see the IFU for more information.

#### Order Number Key



[www.cookmedical.com](http://www.cookmedical.com)

**COOK MEDICAL INCORPORATED**  
P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.  
Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

**COOK (CANADA) INC.**  
111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA  
Phone: 905.640.7110, Toll Free: 800.668.0300

**WILLIAM A. COOK AUSTRALIA PTY. LTD.**  
95 Brandl Street, Brisbane Technology Park, Eight Mile Plains  
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

**WILLIAM COOK EUROPE ApS**  
Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2009 PI-DM-ADV14LP-0907-487-01EN

AORTIC INTERVENTION

CARDIOLOGY

CRITICAL CARE

ENDOSCOPY

INTERVENTIONAL RADIOLOGY

PERIPHERAL INTERVENTION

SURGERY

UROLOGY

WOMEN'S HEALTH