




# Advance<sup>®</sup> 35LP

LOW PROFILE PTA BALLOON DILATATION CATHETER



Expand your procedural options with the first  
8 mm x 8 cm balloon able to pass through a 5 Fr sheath.

#### PIONEERING REWRAP TECHNOLOGY

Starcast™ shoulders provide superior sheath compatibility and sheath reinsertion abilities, allowing repeated treatment during the same procedure.

#### TEXTURED BALLOON SURFACE

Pliiform™ technology reduces surface friction; combined with the lowest crossing profile, it eases crossing of the tightest lesions

#### SUPERIOR SHAFT TECHNOLOGY

Cook's optimized D-shaped lumen shaft design ensures rapid deflation while maintaining superb pushability and trackability.




# Advance® 35LP

## LOW PROFILE PTA BALLOON DILATATION CATHETER

Global Product Number	Order Number	Compatible Sheath Fr	Nominal Pressure atm	Rated Burst Pressure atm
<b>80 cm Catheter Length</b>				
G52337	PTA5-35-80-5-2.0	5	10	15
G52338	PTA5-35-80-5-4.0	5	10	15
G52339	PTA5-35-80-5-6.0	5	10	15
G52340	PTA5-35-80-5-8.0	5	10	15
G52341	PTA5-35-80-5-10.0	5	10	15
G52344	PTA5-35-80-6-2.0	5	10	15
G52345	PTA5-35-80-6-4.0	5	10	15
G52346	PTA5-35-80-6-6.0	5	10	15
G52347	PTA5-35-80-6-8.0	5	10	15
G52348	PTA5-35-80-6-10.0	5	10	15
G52350	PTA5-35-80-6-14.0	5	10	15
G52351	PTA5-35-80-7-2.0	5	8	12
G52352	PTA5-35-80-7-4.0	5	8	12
G52353	PTA5-35-80-7-6.0	5	8	12
G52354	PTA5-35-80-7-8.0	5	8	12
G52355	PTA5-35-80-7-10.0	6	8	12
G52356	PTA5-35-80-7-12.0	6	8	12
G52357	PTA5-35-80-7-14.0	6	8	12
G52358	PTA5-35-80-8-2.0	5	8	11
G52359	PTA5-35-80-8-4.0	5	8	11
G52360	PTA5-35-80-8-6.0	5	8	11
G52361	PTA5-35-80-8-8.0	5	8	11
G52362	PTA5-35-80-8-10.0	6	8	11
G52363	PTA5-35-80-8-12.0	6	8	11
G52364	PTA5-35-80-8-14.0	6	8	11
G52365	PTA5-35-80-9-3.0	6	8	11
G52366	PTA5-35-80-9-4.0	6	8	11
G52367	PTA5-35-80-9-6.0	6	8	11
G52368	PTA5-35-80-9-8.0	6	8	11
G52369	PTA5-35-80-9-10.0	6	8	11
G52370	PTA5-35-80-9-12.0	6	8	11
G52371	PTA5-35-80-9-14.0	6	8	11
G52372	PTA5-35-80-10-3.0	6	8	11
G52373	PTA5-35-80-10-4.0	6	8	11
G52374	PTA5-35-80-10-6.0	6	8	11
G52375	PTA5-35-80-10-8.0	6	8	11
G52376	PTA5-35-80-10-10.0	7	8	11
G52377	PTA5-35-80-10-12.0	7	8	11
G52378	PTA5-35-80-10-14.0	7	8	11
G52379	PTA5-35-80-12-4.0	7	8	11
G52380	PTA5-35-80-12-8.0	7	8	11

Global Product Number	Order Number	Compatible Sheath Fr	Nominal Pressure atm	Rated Burst Pressure atm
<b>135 cm Catheter Length</b>				
G52245	PTA5-35-135-3-2.0	5	10	15
G52246	PTA5-35-135-3-4.0	5	10	15
G52247	PTA5-35-135-3-8.0	5	10	15
G52248	PTA5-35-135-3-10.0	5	10	15
G52249	PTA5-35-135-3-12.0	5	10	15
G52250	PTA5-35-135-3-14.0	5	10	15
G52252	PTA5-35-135-4-4.0	5	10	15
G52253	PTA5-35-135-4-8.0	5	10	15
G52255	PTA5-35-135-4-12.0	5	10	15
G52258	PTA5-35-135-5-4.0	5	10	15
G52260	PTA5-35-135-5-8.0	5	10	15
G52262	PTA5-35-135-5-12.0	5	10	15
G52264	PTA5-35-135-6-2.0	5	10	15
G52265	PTA5-35-135-6-4.0	5	10	15
G52266	PTA5-35-135-6-6.0	5	10	15
G52267	PTA5-35-135-6-8.0	5	10	15
G52269	PTA5-35-135-6-12.0	5	10	15
G52270	PTA5-35-135-6-14.0	5	10	15
G52272	PTA5-35-135-7-4.0	5	8	12
G52273	PTA5-35-135-7-6.0	5	8	12
G52274	PTA5-35-135-7-8.0	5	8	12
G52276	PTA5-35-135-7-12.0	6	8	12
G52278	PTA5-35-135-8-2.0	5	8	11
G52279	PTA5-35-135-8-4.0	5	8	11
G52280	PTA5-35-135-8-6.0	5	8	11
G52281	PTA5-35-135-8-8.0	5	8	11
G52282	PTA5-35-135-8-10.0	6	8	11
G52283	PTA5-35-135-8-12.0	6	8	11
G52284	PTA5-35-135-8-14.0	6	8	11
G52286	PTA5-35-135-9-4.0	6	8	11
G52287	PTA5-35-135-9-6.0	6	8	11
G52288	PTA5-35-135-9-8.0	6	8	11
G52289	PTA5-35-135-9-10.0	6	8	11
G52293	PTA5-35-135-10-4.0	6	8	11
G52294	PTA5-35-135-10-6.0	6	8	11
G52295	PTA5-35-135-10-8.0	6	8	11
G52297	PTA5-35-135-10-12.0	7	8	11

### Order Number Key



The Advance 35LP Low Profile PTA Balloon Dilatation Catheter has been designed for percutaneous transluminal angioplasty (PTA) of lesions in peripheral arteries including iliac, renal, popliteal, infrapopliteal, femoral and iliofemoral as well as obstructive lesions of native or synthetic arterio-venous dialysis fistulae.

Please see IFU for more information.



[www.cookmedical.com](http://www.cookmedical.com)

### COOK MEDICAL INCORPORATED

P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.  
Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

### COOK (CANADA) INC.

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA  
Phone: 905.640.7110, Toll Free: 800.668.0300

### WILLIAM A. COOK AUSTRALIA PTY. LTD.

95 Brandl Street, Brisbane Technology Park, Eight Mile Plains  
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

### WILLIAM COOK EUROPE ApS

Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2009

PIL-2M-ADV35LPDS-EU-200912

AORTIC INTERVENTION

CRITICAL CARE

ENDOSCOPY

INTERVENTIONAL RADIOLOGY

LEAD MANAGEMENT

PERIPHERAL INTERVENTION

SURGERY

UROLOGY

WOMEN'S HEALTH