

# Measurement / Device Selection Form

## Patient Confidential Information

The following information is required to ensure that the appropriate devices and backups are available for the procedure.


Patient ID:

Physician:

Institution:

Imaging Date:

## Measurements – Diameter measurements are inner wall to inner wall



## NOTES:

# Device Selection Patient Confidential Information

The following information is required to ensure that the appropriate devices and backups are available for the procedure.

Patient ID:

Institution:

Physician:

Imaging Date:

Intended Trunk-Ipsilateral Leg Prosthesis Introduction Site:

Right  Left

NOTE: Recommended prosthesis oversizing relative to the aorta is approximately 10–21%. Recommended prosthesis oversizing relative to the iliac artery is approximately 7–25%. All sizes are nominal.

Intended Trunk-Ipsilateral Leg Component: aortic end diameter (mm) x ipsilateral leg diameter (mm) x length (cm)

DEVICE SIZE	QTY.	GORE® C3® DELIVERY SYSTEM*	DEVICE SIZE	QTY.	GORE® C3® DELIVERY SYSTEM*	DEVICE SIZE	QTY.	GORE® C3® DELIVERY SYSTEM*
23 x 12 x 12		RLT231212	26 x 12 x 12		RLT261212	28.5 x 12 x 12		RLT281212
23 x 12 x 14		RLT231214	26 x 12 x 14		RLT261214	28.5 x 12 x 14		RLT281214
23 x 12 x 16		RLT231216	26 x 12 x 16		RLT261216	28.5 x 12 x 16		RLT281216
23 x 12 x 18		RLT231218	26 x 12 x 18		RLT261218	28.5 x 12 x 18		RLT281218
23 x 14.5 x 12		RLT231412	26 x 14.5 x 12		RLT261412	28.5 x 14.5 x 12		RLT281412
23 x 14.5 x 14		RLT231414	26 x 14.5 x 14		RLT261414	28.5 x 14.5 x 14		RLT281414
23 x 14.5 x 16		RLT231416	26 x 14.5 x 16		RLT261416	28.5 x 14.5 x 16		RLT281416
23 x 14.5 x 18		RLT231418	26 x 14.5 x 18		RLT261418	28.5 x 14.5 x 18		RLT281418
						31 x 14.5 x 13		RLT311413
						31 x 14.5 x 15		RLT311415
						31 x 14.5 x 17		RLT311417
						35 x 14.5 x 14		RLT351414
						35 x 14.5 x 16		RLT351416
						35 x 14.5 x 18		RLT351418

Trunk-Ipsilateral Leg Endoprosthesis Sizing Guide:

INTENDED AORTIC INNER DIAMETER (MM)	AORTIC ENDOPROSTHESIS DIAMETER (MM)	INTENDED ILIAC INNER DIAMETER (MM)	ILIAC ENDOPROSTHESIS DIAMETER (MM)	OVERALL DEVICE LENGTH (CM)
19–21	23	10–11 / 12–13.5	12 / 14.5	12 / 14 / 16 / 18
22–23	26	10–11 / 12–13.5	12 / 14.5	12 / 14 / 16 / 18
24–26	28.5	10–11 / 12–13.5	12 / 14.5	12 / 14 / 16 / 18
27–29	31	12–13.5	14.5	13 / 15 / 17
30–32	35	12–13.5	14.5	14 / 16 / 18

\* Products listed may not be available in all markets.

Contralateral Leg Endoprosthesis Sizing Guide:

INTENDED ILIAC INNER DIAMETER (MM)	ILIAC ENDOPROSTHESIS DIAMETER (MM)	CONTRALATERAL LEG LENGTH (CM)	OVERALL DEVICE LENGTH (CM) WHEN INSERTED IN 23, 26, 28.5 MM TRUNKS	OVERALL DEVICE LENGTH (CM) WHEN INSERTED IN 31 MM TRUNKS	OVERALL DEVICE LENGTH (CM) WHEN INSERTED IN 35 MM TRUNKS
10–11	12	10 / 12 / 14	14 / 16 / 18	15 / 17 / 19	16 / 18 / 20
12–13.5	14.5	10 / 12 / 14	14 / 16 / 18	15 / 17 / 19	16 / 18 / 20
13.5–14.5	16**	9.5 / 11.5 / 13.5	13.5 / 15.5 / 17.5	14.5 / 16.5 / 18.5	15.5 / 17.5 / 19.5
14.5–16.5	18**	9.5 / 11.5 / 13.5	13.5 / 15.5 / 17.5	14.5 / 16.5 / 18.5	15.5 / 17.5 / 19.5
16.5–18.5	20**	9.5 / 11.5 / 13.5	13.5 / 15.5 / 17.5	14.5 / 16.5 / 18.5	15.5 / 17.5 / 19.5
18.5–21.5	23**	10 / 12 / 14	14 / 16 / 18	15 / 17 / 19	16 / 18 / 20
21.5–25	27**	10 / 12 / 14	14 / 16 / 18	15 / 17 / 19	16 / 18 / 20

Intended Contralateral Leg Component: iliac end diameter (mm) x length (cm)

QTY.	DEVICE SIZE	CATALOGUE NUMBER	QTY.	DEVICE SIZE	CATALOGUE NUMBER	QTY.	DEVICE SIZE	CATALOGUE NUMBER
	12 x 10	PXC121000		16 x 11.5	PLC161200		20 x 13.5	PLC201400
	12 x 12	PXC121200		16 x 13.5	PLC161400		23 x 10	PLC231000
	12 x 14	PXC121400		18 x 9.5	PLC181000		23 x 12	PLC231200
	14.5 x 10	PXC141000		18 x 11.5	PLC181200		23 x 14	PLC231400
	14.5 x 12	PXC141200		18 x 13.5	PLC181400		27 x 10	PLC271000
	14.5 x 14	PXC141400		20 x 9.5	PLC201000		27 x 12	PLC271200
	16 x 9.5	PLC161000		20 x 11.5	PLC201200		27 x 14	PLC271400

Intended Aortic Extender: for 23, 26, 28.5 mm devices, extension length is 0–1.6 cm; for 32 and 36 mm devices, extension length is 0–2.2 cm

QTY.	INTENDED AORTIC INNER DIAMETER (MM)	ENDOPROSTHESIS DIAMETER (MM)	ENDOPROSTHESIS LENGTH (CM)	CATALOGUE NUMBER
	19–21	23	3.3	PLA230300
	22–23	26	3.3	PLA260300
	24–26	28.5	3.3	PLA280300
	27–29	32	4.5	PLA320400
	30–32	36	4.5	PLA360400

Intended Iliac Extender:\*\* extension length range: 0–4 cm, total extender length: 7 cm

QTY.	INTENDED ILIAC INNER DIAMETER (MM)	ENDOPROSTHESIS DIAMETER (MM)	CATALOGUE NUMBER
	8–9	10	PXL161007
	10–11	12	PXL161207
	12–13.5	14.5	PXL161407

\*\* PLEASE NOTE: All 16 / 18 / 20 / 23 / 27 mm diameter Contralateral Leg Endoprostheses can be used as Iliac Extenders.

Recommended Accessory Selection

GORE® DrySeal Sheath with Hydrophilic Coating GORE® Tri-Lobe Balloon Catheter:

QTY.	DESCRIPTION	CATALOGUE NUMBER	QTY.	INNER VESSEL DIAMETER (MM)	CATALOGUE NUMBER
	12 Fr	DSL1228		16–34	BCM1634
	14 Fr	DSL1428			
	16 Fr	DSL1628			
	18 Fr	DSL1828			

Case planning / procedural notes / comments: